



TECHNICAL DATA SHEET

BOVINE COLLECTION DUMMY WITH HYDRAULIC HEIGHT ADJUSTMENT

Product description

Cow dummy for safe and easy collection of bull semen, with hydraulic height adjustment, mobile or stationary

Properties

- Body made of rigid, zinc-coated steel
- High quality padding
- Hot-dip galvanised body with washable fabric-reinforced PVC skin
- Inside seat for collection technician
- Hydraulic height adjustment for bulls of different sizes
- Comfortable semen collection for the bull
- Forward and backward movements to stimulate the bull (mobile only)
- Separated, smooth hindquarters reduce risk of injury to bull's penis
- Hygienically designed hindquarters for disposable sanitary covers

Dimensions

Dummy measurements:

- Length: 2506 mm
- Width: 974 mm (stationary), 1293 mm (mobile)
- Height: 1240–1443 mm (mounting side)

Body measurements: approx. 2200 x 595 mm (L x W)



Order information

Bovine collection dummy, mobile,
rubber tires with security breaks

REF. : [11000/1000](#)

Bovine collection dummy, stationary

REF. : [11000/0500](#)

Accessories

Disposable sanitary cover, 100/package

REF. : [11000/0200](#)

Plastic cover for dummy cow,
with leather belt

REF. : [11000/0100](#)