

ADAPT HOLDING with HA



ABT's ADAPT HOLDING with HA medium provides embryos a sustainable and healthy environment after collection and before transfer or freezing. We've adapted our ABT Holding formulation by excluding components of animal origin, helping to alleviate concerns about temperature controlled shipping and storage as well as potential bio-hazards. ADAPT Holding contains hyaluronan (HA), a macro-molecule which replaces Bovine Serum Albumin commonly used in other embryo transfer media.

While ADAPT Holding is not meant to act as an extended culture medium for incubation, it does contain a variety of components such as amino acids, nutrients and energy sources that facilitate embryo viability. Embryos can be transferred using a straw containing ADAPT Holding with HA. For specific seeding instructions please read the product insert included in each shipment.

ADAPT Holding starts with ultra-purified Type 1 water. All ABT products are manufactured in a cleanroom environment and go through strict quality control practices, ensuring consistency and quality from lot to lot.

- Includes essential and non-essential amino acids, energy sources, Gentamycin as an antibiotic and is Hepes buffered
- 100% traceable components
- Filtered to 0.1µm; no additional filtering required
- Tight quality control standards equal consistent results
- Available in durable EVA bags and 7ml LDPE Single Use Dropper Bottle
- Made in the USA with more than two decades of experience

PRODUCT #	SIZE
ADA024	50 ml
ADA724	7 ml, Pack of 5

Storage: 2 - 30° C
Shelf Life: 18 months from date of manufacture

ABOUT ABT 360

Located in Pullman, Washington, USA, ABT 360 is owned and operated by Kathy Shearer who has 25 years of experience manufacturing embryo transfer media and Alyssa Hood who has assisted in and managed that same manufacturing facility since 1995. Formed with integrity and experience, ABT 360 has what it takes to support all of your embryo transfer media needs. We look forward to continuing our relationships with the veterinary reproduction community and we will strive to provide quality and consistency in our products, supported by exceptional customer service.



1615 NE EASTGATE BLVD, SEC H
PULLMAN, WA 99163 USA
Phone: 509-592-8144
www.abt360llc.com
alyssa@abt360llc.com
facebook.com/ABT360
instagram.com/ABT_360