



Bovine artificial insemination

QuickLock AI gun

- All-purpose gun for all types of straws and sheaths, split or non-split
- Stainless steel and sterilisable parts only
- One hand only sheath fixation and release
- Unique design for safe and centered straw alignment every time - no waste of semen
- Separate straw release for control of information after AI
- No loose parts - safety catch prevents the plunger from sliding out
- Identification with different colored stylets
- All parts can be replaced

QuickLock Classic

Sheath fixation and release very easy by just pressing a knob

QuickLock Classic [REF.: 17025/0000](#)

QuickLock 2000

- Ergonomic design
- Only 1 cm of the sheath length is required for fixation

QuickLock 2000 [REF.: 17026/0000](#)


QuickLock Ergo

- Ideal for older cows due to increased gun length
- Smallest ergonomic handle available
- Easiest gun as far as cleaning is concerned

QuickLock Ergo [REF.: 17024/0000](#)

QuickLock stylets

Available in different colors.







Universal sheaths

- Suitable for all types of semen straws
- Absolutely tight fit with all types of straws
- Gentle insemination due to a soft shape of the tip
- Transparent adapters for better visual control of straw position
- 50/bag, 1000/box

Universal sheath split (transparent-green) REF: [17007/0001](#)

Universal sheath non-split (transparent-blue) REF: [17007/0000](#)

Hygienic sheath, 80/roll REF: [17028/0100](#)



Thawing unit MT 30/54

- Electronic temperature control
- Central, built-in thawing chamber with lift
- Accuracy: 0.2°C
- Power supply: 12 V DC (car battery)

35°C REF: [17044/0035](#)

37°C REF: [17044/0037](#)

38°C REF: [17044/0038](#)

50°C REF: [17044/0050](#)

Accessories

Transformer for thawing unit 12 V/220 V REF: [17165/1395](#)

Thermometer for thawing unit REF: [17049/0000](#)





QuickLock heater 4.0

- AI and embryo transfer devices are transported at a constant temperature
- Keeps semen/embryo at optimal temperature from thawing to insemination/transfer
- Target temperature inside the bag is freely adjustable
- Power bank as source of electricity for more flexibility and reliability
- Very robust and easy to clean bag material
- Strap allows different ways of carrying the bag

QuickLock heater 4.0 [REF: 17028/0400](#)

QuickLock heater 4.0 with power bank [REF: 17028/0401](#)

Accessories

Power bank for QuickLock heater 4.0 [REF: 17028/0410](#)

USB car lighter charger for power bank [REF: 17028/0420](#)

USB mains charger for power bank [REF: 17028/0430](#)

Inner bag for QuickLock heater 4.0, 100/package [REF: 17028/0201](#)





Insemination kit

Includes: heavy duty case, electronic thawing unit MT 30/54, QuickLock 2000, 100 Universal sheaths, 80 hygienic sheaths for AI gun, ReprJelly lubricant, speculum Varikon, 1 box AI gloves, MiniCutter, tweezers and clamp

Insemination kit [REF.: 17050/0060](#)

Components also available individually:

Insemination case [REF.: 17050/0001](#)

Insemination gloves, 100/box [REF.: 17080/0015](#)

Varikon speculum [REF.: 17064/0380](#)

Straw tweezers, angular shape [REF.: 17060/0000](#)

ReprJelly, 100 ml [REF.: 11907/0100](#)

MiniCutter [REF.: 17062/0010](#)

Clamp [REF.: 17063/0000](#)

Not included in kit:

Transformer for thawing unit [REF.: 17165/1395](#)

